




F0630

6 x 38.5 mL

F0630



Fixed-Angle Rotor, Aluminum

Major applications: Pelleting subcellular organelles, viruses, bacteria, mitochondria, chloroplasts or algae.

Max. rpm	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
26,200*	59,860	454	6 x 38.5 mL 25.3 x 92 mm 1 x 4 in	231 mL

For use in Allegra® 64R, Allegra X-22 Series, and discontinued Allegra 21 Series centrifuges.

No. 361231. F0630 Fixed-Angle Rotor, for 26,200 rpm* operation. Tubes and bottles not included.

Bottles and Tubes

Tube Style/Material	Nominal Volume per Tube (mL)	Part No.	Quantity	Size (mm)	Required Adapters	Tubes per Adapter	g-Force	k Factor	Maximum Speed
Bottles with Cap Assemblies									
Polycarbonate	26.3	355616	3	25 x 89	—	—	59,860	454	26,200
Bottles with Screw Caps									
Polyallomer	30	363073	6	25.3 x 92	—	—	59,860	454	26,200
Polycarbonate	30	363070	6	25.3 x 92	—	—	59,860	454	26,200
Open Top Tubes									
Polyallomer, Thin-wall	34	326825	50	25 x 76	—	—	28,300	—	18,000
Polyallomer, Thick-wall	32	355642	25	25 x 89	—	—	59,860	454	26,200
Ultra-Clear™	38.5	344058	50	25 x 89	—	—	59,860	454	26,200

* In the Allegra 64R centrifuge the maximum rpm is 15,300 and the maximum g-force is 20,450. In the Allegra X22, the maximum rpm is 14,500 and the maximum g-force is 18,300. In the Allegra X-22R, This rotor is rated at a lower speed in the Allegra 21 and 21R centrifuges, so g-forces and k-factors are different. See the rotor manual for specific information.

Rotor Replacement Parts

369353	Rotor Lid
961923	Rotor Lid O-ring
361367	Tie-down Screw
361371	T-handle Rotor Wrench