




F0650

6 x 50 mL

F0650



Fixed-Angle Rotor, Aluminum

Major applications: Density gradient separations of erythrocytes, cell lysate fractions and granules, as well as different separations of DNA, proteins, and viruses.

Max. rpm	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
21,000	41,400	795	6 x 50 mL 29 x 107 mm 1.125 x 4.25 in	300 mL

For use in the Allegra® 64R centrifuge.

No. 364610. F0650 Fixed-Angle Rotor, for 21,000 rpm operation. Tubes and bottles not included.

Bottles and Tubes

Tube Style/Material	Nominal Volume per Tube (mL)	Part No.	Quantity	Size (mm)	Required Adapters	Tubes per Adapter	g-Force	k Factor	Maximum Speed
Bottles with Cap Assemblies									
Polyallomer	50	357001	6	29 x 104	—	—	41,400	795	21,000
	50	361694	24	29 x 104	—	—	41,400	795	21,000
Polycarbonate	50	357000	6	29 x 104	—	—	41,400	795	21,000
Bottles with Screw Caps									
Polyallomer	50	357003	25	29 x 104	—	—	41,400	795	21,000
Polycarbonate	50	357002	25	29 x 104	—	—	41,400	795	21,000
Teflon* with High-speed Screw Cap	50	363076	8	28.5 x 107	—	—	41,400	795	21,000
Tubes									
Polycarbonate	50	363647	25	29 x 104	—	—	41,400	795	21,000
Polycarbonate, Graduated	50	363075	8	29 x 104	—	—	32,300	—	18,500
Polypropylene	50	357007	25	29 x 104	—	—	41,400	795	21,000
Tubes with Snap-On Caps									
Polycarbonate	50	363664	25	29 x 104	—	—	28,611	973	16,500
Polypropylene	50	357005	25	29 x 104	—	—	28,611	973	16,500

* Teflon is a registered trademark of E. I. Du Pont de Nemours & Company.

Rotor Replacement Parts

369355	Rotor Lid
961930	Rotor Lid O-ring
361367	Tie-down Screw
361371	T-handle Rotor Wrench