




F1010

10 x 10 mL

F1010



Fixed-Angle Rotor, Aluminum

Major applications: Pelleting subcellular organelles, viruses, bacteria, mitochondria, chloroplasts or algae.

Max. rpm	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
26,000*	57,440	500	10 x 10 mL 16.1 x 81.1 mm 0.625 x 3 in	100 mL

For use in Allegra® 64R, Allegra X-22 Series and discontinued Allegra 21 Series centrifuges.

1

No. 361221. F1010 Fixed-Angle Rotor, for 26,000 rpm* operation. Tubes and bottles not included.

Bottles and Tubes

Tube Style/Material	Nominal Volume per Tube (mL)	Part No.	Quantity	Size (mm)	Required Adapters	Tubes per Adapter	g-Force*	k Factor	Maximum Speed*
Bottles with Screw Caps									
Polyallomer	10	364695	10	16.1 x 81.1	—	—	57,400	500	26,000
Polycarbonate	10	355672	25	16 x 80	—	—	57,400	500	26,000
Teflon†	10	364693	10	16.1 x 81.1	—	—	57,400	500	26,000
Tubes									
Polyallomer, Thick-wall	6.5	355646‡	25	16 x 64	—	—	—	—	26,000
	10	355640	25	16 x 76	—	—	57,400	500	26,000
Polycarbonate	6.5	355647‡	25	16 x 64	—	—	—	—	26,000
	10	355630	25	16 x 76	—	—	57,400	500	26,000
Ultra-Clear™	13.5	344085	50	16 x 76	—	—	57,400	500	26,000

* In the Allegra 64R centrifuge. This rotor is rated at 15,300 rpm/19,926 x g in the Allegra X-22R and at 14,500 rpm/17,896 x g in the Allegra X-22.

† Teflon is a registered trademark of E. I. Du Pont de Nemours & Company.

‡ To facilitate tube removal, place rubber pad (part number 342602) at the bottom of the tube cavity before inserting the tube.

Rotor Replacement Parts

369351	Rotor Lid
961923	Rotor Lid O-ring
361367	Tie-down Screw
361371	T-handle Rotor Wrench