


JA-21

18 x 10 mL

JA-21



Fixed-Angle Rotor, Aluminum

Major applications: Rapid, efficient separation of many samples in small volumes, such as viruses, bacteriophage, mitochondria, and nuclei.

Max. RPM	Max. <i>g</i>	<i>k</i> Factor	Number of Tubes Volume/Size	Rotor Capacity
21 000	50 400	470	18 x 10 mL 16 x 80 mm 0.625 x 3 in	180 mL

For use in Avanti® J Series and J2 Series centrifuges.

No. 334845. JA-21 Fixed-Angle Rotor for 21 000 rpm operation.

Tubes and Bottles

Tube Style/Material	Nominal Volume per Tube (mL)	Part No.	Quantity	Size (mm)	Required Adapters	Tubes per Adapter	<i>g</i> -Force	<i>k</i> Factor	Maximum Speed
Bottles with Screw Caps									
Polycarbonate	10.0	355672	25	16 x 80	—	—	50 400	469	21 000
Open-Top Tubes									
Thickwall Polyallomer	10.0	355640	25	16 x 76	—	—	50 400	470	21 000
Thickwall Polycarbonate	10.0	355630	25	16 x 76	—	—	50 400	470	21 000

See chart on page 2–6 for adapters used with non-Beckman Coulter tubes and bottles.

Rotor Replacement Parts

870612	O-ring for Rotor Lid
870980	O-ring for Rotor Hub
343631	Lid Attachment Bolt Assembly
369349	Rotor Lid Assembly (Tie-down Assembly included)