


JA-25.50



8 x 50 mL

JA-25.50



Fixed-Angle Rotor, Aluminum

Major applications: Harvesting bacteria, processing tissue homogenates, subcellular particulates, routine pelleting such as precipitates and phase separations.

Max. RPM	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
25 000	75 600	418	8 x 50 mL 29 x 104 mm 1.125 x 4 in	400 mL

For use in Avanti® J Series and J2 Series centrifuges.

No. 363058. JA-25.50[§] Fixed-Angle Rotor with Biosafety Lid.
No. 363055. JA-25.50 Fixed-Angle Rotor with Single Lock Lid.

Tubes and Bottles

Tube Style/Material	Nominal Volume per Tube (mL)	Part No.	Quantity	Size (mm)	Required Adapters	Tubes per Adapter	g-Force	k Factor	Maximum Speed
Bottles with Screw Caps									
Polyallomer	50.0	357003	25	29 x 104	—	—	75 600	418	25 000
Polycarbonate	10.0	355672	25	16 x 80	361703 [†]	1	63 240	350	25 000
	50.0	357002	25	29 x 104	—	—	75 600	418	25 000
Bottles with Cap Assemblies									
Polyallomer	50.0	361694	24	29 x 104	—	—	75 600	418	25 000
Polycarbonate	50.0	361693	24	29 x 104	—	—	75 600	418	25 000
Tubes with Snap-On Caps									
Thickwall Polycarbonate	50.0	363664*	25	29 x 104	—	—	75 600	418	25 000
Thickwall Polypropylene	1.5	356090	500	11 x 38	361705 [†]	3	50 960	182	25 000
	50.0	357005*	25	29 x 104	—	—	75 600	418	25 000
Open-Top Tubes									
Polyallomer	1.5	357448	500	11 x 38	361705 [†]	3	50 960	182	25 000
Thickwall Polycarbonate	50.0	363647	25	29 x 104	—	—	75 600	418	25 000
Thickwall Polypropylene	50.0	357007	25	29 x 104	—	—	75 600	418	25 000

* Order cap separately (pkg. of 25): 356264 (natural), 357359 (red), 357360 (green), 357361 (yellow), or 357362 (blue).

[†] Set of 2.

[§] This rotor was tested to demonstrate containment of microbiological aerosols under normal operating conditions of the associated Beckman Coulter centrifuge, when used and maintained as instructed. Validation of microbiological containment was done at an independent, third-party facility (CAMR, Porton Down, UK or USAMRIID, Ft. Detrick, MD, USA). Improper use or maintenance may affect seal integrity and, thus, containment.

Rotor Replacement Parts

870612	O-ring for Rotor Lid
011757	Small O-ring for Rotor Lid
363023	Lid Assembly, Biosafety (Tie-down Assembly included)
363053	Lid Assembly, Single Lock (Tie-down Assembly included)

Adapters

See chart on page 2–6 for adapters used with non-Beckman Coulter tubes and bottles.

361703 361705

