




MLA-130

10 x 2.0 mL

MLA-130

1



Fixed-Angle Rotor, Titanium

For use in the Optima™ MAX Series Ultracentrifuges.

Major applications: Rapid differential sedimentation (pelleting) of small particles such as subcellular organelles and viruses.

Max. RPM	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
130 000	1 019 000	8.7	10 x 2.0 mL 11 x 32 mm	20 mL

No. 367114. MLA-130 Fixed-Angle Rotor Package. Includes Rotor, one box of 100 each 344625 Polyallomer Quick-Seal® Tubes, 367310 Spacer, 348305 Tube Rack, 975001 O-ring, 927208 Hemostats, 306812 Spinkote™ Lubricant and 335148 Vacuum Grease.

No. 367120. MLA-130 Fixed-Angle Rotor Assembly.

Tubes and Bottles

Tube Style/Material	Vol. (mL)	Part No.	Quantity	Size (mm)	Floating Spacers (qty. 1)	Adapters (qty. 1)	g-Force	k factor	Max. Speed
Quick-Seal®	1.5	344624	100	11 x 25	367310	—	1 008 000	7.0	1 008 000
Polyallomer	2.0	344625	100	11 x 32	367310	—	1 008 000	8.7	1 008 000
Thickwall Polyallomer	1.0	347287	100	11 x 34	—	—	390 000	23	80 000
Thickwall Polycarbonate	1.0	343778	100	11 x 34	—	—	1 008 000	5.8	1 008 000

Rotor Accesories

358312	Quick-Seal Cordless Tube Topper Kit, 60 Hz
358313	Quick-Seal Cordless Tube Topper Kit, 50 Hz (Europe)
358314	Quick-Seal Cordless Tube Topper Kit, 50 Hz (Great Britain)
358315	Quick-Seal Cordless Tube Topper Kit, 50 Hz (Australia)
367803	Quick-Seal Cordless Tube Topper Kit, 50 Hz (Canada)
349387	Tube Topper Rack
361668	Tube Removal Tool
927208	Hemostat, 6 in. curved
961519	Hemostat, 6 in. straight
342025	Fraction Recovery System
347828	Fraction Recovery System Adapter Kit for TL-series tubes
306812	Spinkote Lubricant (2 oz)
335148	Silicone Vacuum Grease
339558	Rotor Cleaning Kit
347404	Rotor Cleaning Brush
339555	Beckman Solution 555 (1 qt)

Rotor Replacement Parts

367116	Rotor Lid
975001	O-ring for Lid
347373	Rotor Vise

Adapters/Spacers

367310

