


MLA-150

8 x 2.0 mL

MLA-150



Fixed-Angle Rotor, Titanium

For use in the Optima™ MAX Series Ultracentrifuges.

Major applications: Rapid differential sedimentation (pelletting) of small particles such as subcellular organelles and viruses.

Max. RPM	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
150 000	1 003 000	10.4*	8 x 2.0 mL 11 x 32 mm	16 mL

* Based on cell hole dimensions. k-factor = 4.58 if run with 1 mL sample and open top tube.

No. 393490. MLA-150 Rotor Package. Includes Rotor, one box of 100 each 343778 Polycarbonate tubes, one box of 50 each 344625 Quick-Seal® Polyallomer tubes, 348305 Tube Rack, 975001 O-ring, 927208 Hemostats, 306812 Spinkote™ Lubricant and 335148 Vacuum Grease.

No. 393489. MLA-150 Fixed-Angle Rotor Assembly.

Tubes and Bottles

Tube Style/Material	Vol. (mL)	Part No.	Quantity	Size (mm)	Floating Spacers (pkg. of 8)	Adapters (qty. 1)	g-Force	k factor	Max. Speed
Quick-Seal® Polyallomer	1.5	344624	50	11 x 25	393570	—	989 000	6.2	150 000
Polyallomer	2.0	344625	50	11 x 32	393570	—	989 000	8.0	150 000
Thickwall Polyallomer	1.0	347287	100	11 x 34	—	—	387 000	11.4	95 000
Thickwall Polycarbonate	1.0	343778	100	11 x 34	—	—	964 000	4.6	150 000

Rotor Accessories

358312 Quick-Seal Cordless Tube Topper Kit, 60 Hz
 358313 Quick-Seal Cordless Tube Topper Kit, 50 Hz (Europe)
 358314 Quick-Seal Cordless Tube Topper Kit, 50 Hz (Great Britain)
 358315 Quick-Seal Cordless Tube Topper Kit, 50 Hz (Australia)
 367803 Quick-Seal Cordless Tube Topper Kit, 50 Hz (Canada)
 349387 Tube Topper Rack
 361668 Tube Removal Tool
 927208 Hemostat, 6 in. curved
 961519 Hemostat, 6 in. straight
 342025 Fraction Recovery System
 347828 Fraction Recovery System Adapter Kit for TL-series tubes
 306812 Spinkote Lubricant (2 oz)
 335148 Silicone Vacuum Grease
 339558 Rotor Cleaning Kit
 347404 Rotor Cleaning Brush
 339555 Beckman Solution 555 (1 qt)

Rotor Replacement Parts

860218 O-ring for Lid
 393485 Rotor Lid Assembly

Adapters/Spacers

344674

