


MLA-80

8 x 8 mL

MLA-80



Fixed-Angle Rotor, Titanium

For use in Optima™ MAX Series high-capacity personal micro-ultracentrifuges.

Major applications: Large volume differential sedimentation (pelletting) of subcellular organelles and viruses.

Max. RPM	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
80 000	444 000	29	8 x 8 mL 16 x 64 mm	64 mL

No. 367096. MLA-80 Fixed-Angle Rotor Package. Includes Rotor, one box of 50 each 344621 Tubes, 367094 Quick-Seal® Bell Top Spacer, 892292 O-ring, 348123 Tube Rack, 927208 Hemostats, 338765 Stem Lifter, 306812 Spinkote™ Lubricant and 335148 Vacuum Grease.

No. 367087. MLA-80 Fixed-Angle Rotor Assembly.

Tubes and Bottles

Tube Style/Material	Vol. (mL)	Part No.	Quantity	Size (mm)	Required Caps/Spacers (qty. 1)	Required Adapters (qty. 1)	g-Force	k factor	Max. Speed
Quick-Seal®	4.2	356562	50	16 x 38	367094	—	444 000	18	80 000
Polyallomer	6.3	345830	50	16 x 45	367094	—	444 000	23	80 000
	8.0	344621	50	16 x 58	367094	—	444 000	29	80 000
Thickwall Polyallomer	6.5	355646	25	16 x 64		—	84 900	153	35 000
Thickwall Polycarbonate	6.5	355647	25	16 x 64		—	417 000	19	80 000

Rotor Accessories

358312 Quick-Seal Cordless Tube Topper Kit, 60 Hz
 358313 Quick-Seal Cordless Tube Topper Kit, 50 Hz (Europe)
 358314 Quick-Seal Cordless Tube Topper Kit, 50 Hz (Great Britain)
 358315 Quick-Seal Cordless Tube Topper Kit, 50 Hz (Australia)
 367803 Quick-Seal Cordless Tube Topper Kit, 50 Hz (Canada)
 348123 Tube Topper Rack
 361668 Tube Removal Tool
 338765 Floating Spacer Removal Tool
 927208 Hemostat, 6 in. curved
 961519 Hemostat, 6 in. straight
 306812 Spinkote Lubricant (2 oz)
 335148 Silicone Vacuum Grease
 339558 Rotor Cleaning Kit
 347404 Rotor Cleaning Brush
 339555 Beckman Solution 555 (1 qt)

Rotor Replacement Parts

367089 Rotor Lid
 892292 O-ring
 346133 Rotor Vise

Adapters/Spacers

367094

