


SW 28

6 x 38.5 mL

SW 28



Swinging-Bucket Rotor, Aluminum Head, and Titanium Buckets

For use in instruments classified: **H R S**

Major applications: Differential centrifugation of subcellular fractions and viruses.

Max. RPM	Max. g	k Factor	Number of Tubes Volume/Size	Rotor Capacity
28 000	141 000	246	6 x 38.5 mL 1 x 3½ in 25 x 89 mm	231 mL

No. 342204. SW 28 Rotor Package with 2 boxes 326823 Polyallomer Tubes, 2 boxes 344058 Ultra-Clear™ Tubes, 331186 Bucket Holder Rack, 306812 Spinkote™ Lubricant, 335148 Vacuum Grease, 342211 Overspeed Disk, and 812715 O-rings.

No. 342207. SW 28 Rotor Assembly.

Tubes									
Tube Style/Material	Vol. (mL) <i>g-Max</i>	Part No.	Quantity	Size (mm)	Floating Spacers (qty. 1)	Adapters (qty. 1)	g-Force	k Factor	Max. Speed
OptiSeal Polyallomer	32.4	361625	56	25 x 77	361669*	—	141 000	206	28 000
Quick-Seal®	15.0	343664	50	25 x 38	355536	—	141 000	87	28 000
Polyallomer	27.0	343665	50	25 x 64	355536	—	141 000	164	28 000
Bell-Top	33.5	344623	50	25 x 83	355536	—	141 000	234	28 000
Quick-Seal	8.4	358652	50	25 x 38	355536	358156**	141 000	87	28 000
Konical™	22.5	358654	50	25 x 76	355536	358156**	141 000	206	28 000
Polyallomer	28.0	358651	50	25 x 83	355536	358156**	141 000	234	28 000
Bell-Top									
Konical	25.5	358125	50	25 x 76	—	358156**	141 000	206	28 000
Polyallomer	31.5	358126	50	25 x 89	—	358156**	141 000	246	28 000
Thickwall	31.0	355642	25	25 x 89	—	—	141 000	246	28 000
Polyallomer									
Thinwall	38.5	326823	50	25 x 89	—	—	175 000	204	32 000
Polyallomer									
Thickwall	31.0	355631	25	25 x 89	—	—	141 000	246	28 000
Polycarbonate									
Ultra-Clear™	38.5	344058	50	25 x 89	—	—	141 000	246	28 000

* Set of 2.

** Package of 6.

Note: To help you locate g-Max tubes in the rotor listings, these tubes are highlighted with bold, italic type.

Rotor Replacement Parts

331186 Bucket Holder Rack
332400 Rotor Stand
342211 Overspeed Disk, 28 000 rpm
342217 Bucket Assembly, Titanium, with Caps and O-rings, matched set of 6
812715 O-ring for Bucket
342179 Bucket Cap (each)

Extra Buckets for the SW 28

The long slender buckets of the SW 28.1 can be used with the SW 28 rotor head.

342212 Bucket Assembly SW 28.1, Titanium, with Caps and O-rings, matched set of 6

Adapters/Spacers

355536

358156

361669

