




TA-10-250

6 x 250 mL

TA-10-250



Fixed-Angle Rotor, Aluminum

Major applications: Pelleting cells from large volumes or cell particles from tissue homogenates. Short-column methods (such as partially filled tubes) may also be used to purify large quantities of virus in a cushion gradient.

Max. rpm	Max. g	Rotor Capacity	Approximate Accel/Decel Time (min:sec)
10,000	15,300	6 x 250 mL 62 x 141 mm 2.5 x 5.75 in.	2:00/2:00

For use in the Allegra® 25R centrifuge.

No. 368293. TA-10-250 Fixed-Angle Rotor.

Bottles and Tubes

Tube Style/Material	Nominal Volume per Tube (mL)	Part No.	Quantity	Size (mm)	Required Adapters	Tubes per Adapter	g-Force	k Factor	Maximum Speed
Bottles with Cap Assemblies									
Polyallomer	50	357001	6	29 x 104	356997	1	—	—	10,000
	50	361694	24	29 x 104	356997	1	—	—	10,000
Polypropylene	50	357000	6	29 x 104	356997	1	—	—	10,000
Polypropylene, Wide-mouth	250	356011	6	62 x 122	—	—	15,300	3,450	10,000
Polycarbonate, Wide-mouth	250	356013	6	62 x 122	—	—	15,300	3,450	10,000
Bottles with Screw Caps									
Polyallomer	50	357003	25	29 x 104	356997	1	—	—	10,000
Polycarbonate	50	357002	25	29 x 104	356997	1	—	—	10,000
Tubes with Snap-On Caps									
Polycarbonate	50	363664	25	29 x 104	356997	1	—	—	10,000
Polypropylene	50	357005	25	29 x 104	356997	1	—	—	10,000
Conical Tubes									
Polycarbonate with Cap	230	356987	6	62 x 141	356983	1	—	—	10,000
Polypropylene with Cap	230	356989	6	62 x 141	356983	1	—	—	10,000
Open-Top Tubes									
Polycarbonate	50	363647	25	29 x 104	356997	1	—	—	10,000
Polypropylene	50	357007	25	29 x 103	356997	1	—	—	10,000

Rotor Replacement Parts

368245	Tie-down Screw
368246	T-handle Rotor Wrench

Adapters

356983

356997

