




TS-5.1-500

Biosafety

4 x 500 mL

TS-5.1-500

1



Swinging Bucket Rotor, Aluminum

Major applications: Sedimenting protein precipitates, large particles, cells and cell debris.

Max. rpm	Max. g	Rotor Capacity	Approximate Accel/Decel Time (min:sec)
5,100	5,540	4 x 500 mL 85 x 135 mm 3.5 x 5.5 in.	0:59/0:47

For use in the Allegra® 25R centrifuge.

No. 368308. TS-5.1-500 Swinging-Bucket Biosafety Rotor.*

Bottles and Tubes

Tube Style/Material	Nominal Volume per Tube (mL)	Part No.	Quantity	Size (mm)	Required Adapters	Tubes per Adapter	g-Force	k Factor	Maximum Speed
Bottles with Cap Assembly									
Polycarbonate	250	356013	6	62 x 122	368456	1	—	—	5,100
Polypropylene	250	356011	6	62 x 122	368456	1	—	—	5,100
Bottles with Screw Caps									
Polyallomer	10	364695	10	16.1 x 81.1	368468	12	—	—	5,100
Polycarbonate	10	355672	25	16 x 80	368468	12	—	—	5,100
	500	368454	2	85 x 135	—	—	5,540	—	5,100
Polypropylene	500	368453	2	85 x 135	—	—	5,540	—	5,100
Tubes with Caps									
Polyallomer	1.5	357448	500	11 x 38	368470	16	—	—	5,100
Open-Top Tubes									
Polyallomer, Thick-wall	10	355640	25	16 x 76	368468	12	—	—	5,100
Polycarbonate	6.5	355647	25	16 x 64	368468	12	—	—	5,100
	10	355630	25	16 x 76	368468	12	—	—	5,100

* With sealed bucket covers.

Rotor Replacement Parts

368452	Buckets (set of 2)
368472	Bucket Cover (set of 2)
368455	Bucket O-ring (set of 8)
368451	Microplate Carrier (set of 2)
368245	Tie-down Screw
368246	T-handle Wrench

Adapters

